

# bastion



Latex Powder Free  White Gloves



Ultimate solutions in hand protection

# Latex Powder Free White Gloves

<b>Product:</b>	Latex Polymer Coated White Powder Free Gloves
<b>Material:</b>	Natural Rubber Latex
<b>Type:</b>	Non Sterile
<b>Design:</b>	Ambidextrous
<b>Feature:</b>	Fully Textured
<b>Cuff:</b>	Beaded
<b>Usage:</b>	Single Use only
<b>Colour:</b>	White
<b>Packaging:</b>	100 ± 2 pcs/inner, 10 inners/carton

## Dimensions:

Code	Size	Carton Quantity	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG2821	XS	1000	70 - 79	Min. 240	Min. 220
BNG2822	S	1000	80 - 89	Min. 240	Min. 230
BNG2823	M	1000	90 - 99	Min. 240	Min. 230
BNG2824	L	1000	100 - 109	Min. 240	Min. 230
BNG2825	XL	1000	110 - 119	Min. 240	Min. 230

## Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.15
Palm	Min. 0.15
Cuff	Min. 0.10

## Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 18	Min. 14
Elongation (%)	Min. 650	Min. 500

**Powder Level:** Powder free

- Silicone free

