

bastion



Premium Latex  White Powder Free Gloves



Ultimate solutions in hand protection

Premium Latex White Powder Free Gloves

Product:	Premium Latex Powder Free White Gloves
Material:	Natural Rubber Latex
Type:	Non Sterile
Design:	Ambidextrous
Feature:	Fully Textured
Cuff:	Beaded
Usage:	Single Use only
Colour:	White
Packaging:	100 ± 2 pcs/inner, 10 inners/carton

Dimensions:

Code	Size	Carton Quantity	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG2831	XS	1000	70 - 79	Min. 240	Min. 220
BNG2832	S	1000	80 - 89	Min. 240	Min. 230
BNG2833	M	1000	90 - 99	Min. 240	Min. 230
BNG2834	L	1000	100 - 109	Min. 240	Min. 230
BNG2835	XL	1000	110 - 119	Min. 240	Min. 230

Thickness:

Location of Thickness Measurement	Single Wall (mm)		
Finger	Min. 0.15		
Palm	Min. 0.15		
Cuff	Min. 0.10		

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 18	Min. 14
Elongation (%)	Min. 650	Min. 500

Powder Level: Powder Free

- Silicone free



TGA