


bastion



Latex™ High Risk Long Cuff  Blue
Powder Free Gloves



Ultimate solutions in hand protection

Latex™ High Risk Long Cuff Blue Powder Free Gloves

Product:	Latex High Risk Long Cuff Blue Powder Free Gloves
Material:	Natural Rubber Latex
Type:	Non Sterile
Design:	Ambidextrous
Feature:	Fully Textured
Cuff:	Beaded
Usage:	Single Use only
Colour:	Blue
Packaging:	S - XL - 50 ± 2 pcs/inner, 10 inners/carton XXL - 45 ± 2 pcs/inner, 10 inners/carton

Dimensions:

Code	Size	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG2962	S	83 - 88	Min. 300	Min. 295
BNG2963	M	93 - 98	Min. 300	Min. 295
BNG2964	L	103 - 109	Min. 300	Min. 295
BNG2965	XL	115 - 120	Min. 300	Min. 295
BNG2966	XXL	122 - 127	Min. 300	Min. 295

Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.39
Palm	Min. 0.33
Cuff	Min. 0.20

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 18	Min. 14
Elongation (%)	Min. 650	Min. 500
Force at Break (N)	Min. 9	Min. 6

Powder Level:	Powder Free
Compliance Standards:	ASTM D3578-05



TGA