

bastion

Vinyl  Red Powder Free Gloves



Ultimate solutions in hand protection

Vinyl Red Powder Free Gloves

Product:	Vinyl Red Powder Free Gloves
Material:	Polyvinyl Chloride
Type:	Non Sterile
Design:	Ambidextrous
Cuff:	Beaded
Usage:	Single Use only
Colour:	Red
Packaging:	100 pcs/inner, 10 inners/carton S - XL 90 pcs/inner, 10 inners/carton XXL

Dimensions:

Code	Size	Carton Quantity	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG3822	S	1000	80 - 89	Min. 240	Min. 220
BNG3823	M	1000	90 - 99	Min. 240	Min. 230
BNG3824	L	1000	100 - 109	Min. 240	Min. 230
BNG3825	XL	1000	110 - 119	Min. 240	Min. 230
BNG3826	XXL	900	120 - 127	Min. 240	Min. 230

Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.10
Palm	Min. 0.09
Cuff	Min. 0.06

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 11	Min. 11
Elongation (%)	Min. 300	Min. 300

Powder Level: Powder free

- DOP and DEHP free
- Latex Free



TGA