


bastion



Nitrile  **ExtraTough** Green Powder Free
Gloves



Ultimate solutions in hand protection

Nitrile ExtraTough Green Powder Free Gloves

Product: Nitrile ExtraTough Green Powder Free Gloves
Material: Superior Formulation Nitrile
Type: Non Sterile
Design: Ambidextrous
Feature: Micro Textured Fingers
Cuff: Beaded
Usage: Single Use only
Colour: Green
Packaging: S - XL 100/inner, 10 inners/carton
XXL 90/inner, 10 inners/carton

Dimensions:

Code	Size	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG6842	S	70 - 77	Min. 240	Min. 220
BNG6843	M	80 - 87	Min. 240	Min. 230
BNG6844	L	90 - 97	Min. 240	Min. 230
BNG6845	XL	100 - 107	Min. 240	Min. 230
BNG6846	XXL	110 - 117	Min. 240	Min. 230

Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.15
Palm	Min. 0.12
Cuff	Min. 0.08

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400
Force at break (N)	> 6	> 6

Powder Level: Powder Free



TGA