

bastion



Nitrile *UltraTough*®  Black Powder Free Gloves



Ultimate solutions in hand protection

Nitrile *UltraTough*® Black Powder Free Gloves

Product:	Nitrile UltraTough Black Powder Free Gloves
Material:	Nitrile - Latex Free
Type:	Non Sterile
Design:	Ambidextrous
Feature:	Micro Textured
Cuff:	Beaded
Usage:	Single Use only
Colour:	Black
Packaging:	100 pcs/inner, 10 inners/carton (S - XL) 90 pcs/inner, 10 inners/carton (XXL)

Dimensions:

Code	Size	Palm Width (mm)	MGM Length (mm)
BNG6882	S	70 - 79	Min. 240
BNG6883	M	80 - 89	Min. 240
BNG6884	L	90 - 99	Min. 240
BNG6885	XL	100 - 109	Min. 240
BNG6886	XXL	110 - 119	Min. 240

Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.20
Palm	Min. 0.17
Cuff	Min. 0.10

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400



EN455