

# bastion



SensiSafe  Aloe Vera Nitrile Gloves



Ultimate solutions in hand protection

# SensiSafe Aloe Vera Nitrile Gloves

**Product:** SensiSafe Aloe Vera Nitrile Gloves  
**Material:** Nitrile - Sulphur Free, Accelerator Free  
**Type:** Non Sterile  
**Design:** Ambidextrous  
**Feature:** Micro Textured Fingers  
**Cuff:** Beaded  
**Usage:** Single Use only  
**Colour:** Green  
**Packaging:** XS - XL 200/inner, 10 inners/carton

## Dimensions:

Codes	Size	Palm Width (mm)	Length (mm)
BNG7141	XS	70 - 79	Min. 240
BNG7142	S	80 - 89	Min. 240
BNG7143	M	90 - 99	Min. 240
BNG7144	L	100 - 109	Min. 240
BNG7145	XL	110 - 119	Min. 240

## Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.11
Palm	Min. 0.08
Cuff	Min. 0.06

## Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400



EN374  
EN455

ASTM D6319