


bastion



Nitrile Ultra Soft  Black Powder Free Gloves



Ultimate solutions in barrier protection

Nitrile Ultra Soft Black Powder Free Gloves

Product:	Nitrile Ultra Soft Black Powder Free Gloves
Material:	Nitrile
Type:	Non Sterile
Design:	Ambidextrous
Feature:	Micro Textured Fingers
Cuff:	Beaded
Usage:	Single Use only
Colour:	Black
Packaging:	100 pcs/inner, 10 inners/carton

Dimensions:

Code	Size	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG7851	XS	70 - 79	Min. 240	Min. 220
BNG7852	S	80 - 89	Min. 240	Min. 230
BNG7853	M	90 - 99	Min. 240	Min. 230
BNG7854	L	100 - 109	Min. 240	Min. 230
BNG7855	XL	110 - 119	Min. 240	Min. 230

Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.13
Palm	Min. 0.08
Cuff	Min. 0.05

Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400

Powder Level: Powder Free



TGA