

# bastion



Nitrile *UltraSoft*<sup>®</sup>  Blue Powder Free Gloves



Ultimate solutions in hand protection

# Nitrile *UltraSoft*<sup>®</sup> Blue Powder Free Gloves

<b>Product:</b>	Nitrile UltraSoft <sup>®</sup> Blue Powder Free Gloves
<b>Material:</b>	Nitrile
<b>Type:</b>	Non Sterile
<b>Design:</b>	Ambidextrous
<b>Feature:</b>	Micro Textured Fingers
<b>Cuff:</b>	Beaded
<b>Usage:</b>	Single Use only
<b>Colour:</b>	Blue
<b>Packaging:</b>	200 ± 2 pcs/inner, 10 inners/carton

## Dimensions:

Code	Size	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BNG7861	XS	70 - 79	Min. 240	Min. 220
BNG7862	S	80 - 89	Min. 240	Min. 230
BNG7863	M	90 - 99	Min. 240	Min. 230
BNG7864	L	100 - 109	Min. 240	Min. 230
BNG7865	XL	110 - 119	Min. 240	Min. 230
BNG7866	XXL	120 - 129	Min. 240	Min. 230

## Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.13
Palm	Min. 0.08
Cuff	Min. 0.05

## Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400

**Powder Level:** Powder Free



# TGA