



Head Wear - PP Bouffant Caps - 21"



Head Wear - PP Bouffant Caps - 21"

- 14 gsm non woven polypropylene
- Lightweight and comfortable
- For use in operating theatres and for general ward use
- Protects staff and patients from transfer of micro-organisms and particulate matter
- High quality ensuring containment of hair and skin shed
- HACCP certified and TGA listed

PP Bouffant Caps - 21"

Code	Colour	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BNR21124	White	10	100	1000
BNR21134	Blue	10	100	1000
BNR21144	Red	10	100	1000
BNR21154	Green	10	100	1000
BNR21164	Yellow	10	100	1000