

bastion



Gowns



Gowns

- 40 gsm non woven polypropylene
- Excellent protection
- Waist ties
- Velcro closure at neck
- Elastic cuffs
- One size fits all
- TGA listed

PP Clinical Gowns

Code	Colour	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BNR25121	White	100	1	100
BNR25131	Blue	100	1	100

PP Isolation Gowns

Code	Colour	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BNR25261	Yellow	10	10	100